



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesell-
schaft mbH
Gutenbergring 5, 63110 Rodgau



Page 1/15

Safety Data Sheet

According to Annex II to REACH - Regulation (EU) 2020/878 and to Annex II to UK REACH

SECTION 1. Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier			
Product name	OMNI ALGINAT BLAU PFEFFERMINZ		
1.2. Relevant identified uses of the substance or mixture and uses advised against			
Identified Uses	Industrial	Professional	Consumer
Non-chromatic alginate for dental impressions.	-		-
Uses Advised Against			
All uses not indicated in the identified uses.			
1.3. Details of the supplier of the safety data sheet			
Name	LASCOD S.p.A.		
Full address	VIA LUIGI LONGO 18		
District and Country	50019 SESTO FIORENTINO (FI)		
	Italia		
	tel. +39 055 421 5768		
e-mail address of the competent person			
responsible for the Safety Data Sheet	info@lascod.it		
Supplier:	LASCOD S.p.A.		
1.4. Emergency telephone number			
For urgent inquiries refer to			

(AT): +43 14064343 - (BE): +32 070245245 - (BG): +359 29154233 - (CY): 1401



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 2/15

(DK): +45 82121212 - (ES): +34 915620420 - (HR): +385 12348342
(EE): +372 7943794, 16662 - (FI): +358 800 147111 - (FR): +33 (0) 145425959
(DE): 112 - (EL): +30 2107793777 - (IE): +353 (01) 8092166 - (IS): +354 5432222, 112
(LV): +371 67042473, 112 - (LT): +370 (85) 2362052 - (LU): +352 80025500
(MT): 2545 0000 - (NO): +47 22591300 - (NL): +31 (0) 887558000 - (PL): 112
(PT): +351 800250250 - (CZ): +420 224919293 - (RO): +40 212308000
(SK): +421 254774166 - (SI): 112 - (SE): 112 - (HU): +36 80201199

SECTION 2. Hazards identification

2.1. Classification of the substance or mixture

The product is classified as hazardous pursuant to the provisions set forth in (EC) Regulation 1272/2008 (CLP) (and subsequent amendments and supplements). The product thus requires a safety datasheet that complies with the provisions of (EU) Regulation 2020/878.

Any additional information concerning the risks for health and/or the environment are given in sections 11 and 12 of this sheet.

Hazard classification and indication:

Specific target organ toxicity - repeated exposure, category 2	H373	May cause damage to organs through prolonged or repeated exposure.
Eye irritation, category 2	H319	Causes serious eye irritation.

2.2. Label elements

Hazard labelling pursuant to EC Regulation 1272/2008 (CLP) and subsequent amendments and supplements.

Hazard pictograms:						
Signal words:	Warning					
Hazard statements:						
H373	May cause damage to organs through prolonged or repeated exposure.					
H319	Causes serious eye irritation.					
Precautionary statements:						
P280	Wear eye protection / face protection.					
P337+P313	If eye irritation persists: Get medical advice / attention.					



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 3/15

Contains:	DIATOMITE

2.3. Other hazards

On the basis of available data, the product does not contain any PBT or vPvB in percentage $\geq$ than 0,1%.
The product does not contain substances with endocrine disrupting properties in concentration $\geq$ 0.1%.

SECTION 3. Composition/information on ingredients

3.2. Mixtures

Contains:

Identification	x = Conc. %	Classification (EC) 1272/2008 (CLP)	
DIATOMITE			
INDEX -	$70 \leq x < 85$	STOT RE 2 H373	
EC 272-489-0			
CAS 68855-54-9			
REACH Reg. 01-2119488518-22			
DIPOTASSIUM HEXAFLUOROTITANATE			
INDEX -	$0,5 \leq x < 2$	Acute Tox. 4 H302, Eye Dam. 1 H318	
EC 240-969-9		LD50 Oral: 324 mg/kg	
CAS 16919-27-0			
REACH Reg. 01-2119978268-20			

The full wording of hazard (H) phrases is given in section 16 of the sheet.

SECTION 4. First aid measures



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 4/15

4.1. Description of first aid measures

EYES: Remove contact lenses, if present. Wash immediately with plenty of water for at least 15 minutes, opening the eyelids fully. If problem persists, seek medical advice.

SKIN: Remove contaminated clothing. Rinse skin with a shower immediately. Wash contaminated clothing before using it again.

INHALATION: Remove to open air. If the subject stops breathing, administer artificial respiration. Get medical advice/attention immediately.

INGESTION: Get medical advice/attention immediately. Do not induce vomiting. Do not administer anything not explicitly authorised by a doctor.

4.2. Most important symptoms and effects, both acute and delayed

Specific information on symptoms and effects caused by the product are unknown.

4.3. Indication of any immediate medical attention and special treatment needed

Information not available

SECTION 5. Firefighting measures

5.1. Extinguishing media

SUITABLE EXTINGUISHING EQUIPMENT

The extinguishing equipment should be of the conventional kind: carbon dioxide, foam, powder and water spray.

UNSUITABLE EXTINGUISHING EQUIPMENT

None in particular.

5.2. Special hazards arising from the substance or mixture

HAZARDS CAUSED BY EXPOSURE IN THE EVENT OF FIRE

Do not breathe combustion products. The product is combustible and, when the powder is released into the air in sufficient concentrations and in the presence of a source of ignition, it can create explosive mixtures with air. Fires may start or get worse by leakage of the solid product from the container, when it reaches high temperatures or through contact with sources of ignition.

5.3. Advice for firefighters

GENERAL INFORMATION

Use jets of water to cool the containers to prevent product decomposition and the development of substances potentially hazardous for health. Always wear full fire prevention gear. Collect extinguishing water to prevent it from draining into the sewer system. Dispose of contaminated water used for extinction and the remains of the fire according to applicable regulations.

SPECIAL PROTECTIVE EQUIPMENT FOR FIRE-FIGHTERS

Normal fire fighting clothing i.e. fire kit (BS EN 469), gloves (BS EN 659) and boots (HO specification A29 and A30) in combination with self-contained open circuit positive pressure compressed air breathing apparatus (BS EN 137).

SECTION 6. Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesell-
schaft mbH
Gutenbergring 5, 63110 Rodgau



Page 5/15

If there are no contraindications, spray powder with water to prevent the formation of dust.
Wear suitable protective equipment (including personal protective equipment referred to under Section 8 of the safety data sheet) to prevent any contamination of skin, eyes and personal clothing. These indications apply for both processing staff and those involved in emergency procedures.

6.2. Environmental precautions

The product must not penetrate into the sewer system or come into contact with surface water or ground water.

6.3. Methods and material for containment and cleaning up

Collect the leaked product and place it in containers for recovery or disposal. If there are no contraindications, use jets of water to eliminate product residues.

Make sure the leakage site is well aired. Evaluate the compatibility of the container to be used, by checking section 10. Contaminated material should be disposed of in compliance with the provisions set forth in point 13.

6.4. Reference to other sections

Any information on personal protection and disposal is given in sections 8 and 13.

SECTION 7. Handling and storage

7.1. Precautions for safe handling

Before handling the product, consult all the other sections of this material safety data sheet. Avoid leakage of the product into the environment. Do not eat, drink or smoke during use. Remove any contaminated clothes and personal protective equipment before entering places in which people eat.

7.2. Conditions for safe storage, including any incompatibilities

Store only in the original container. Store the containers sealed, in a well ventilated place, away from direct sunlight. Keep containers away from any incompatible materials, see section 10 for details.

7.3. Specific end use(s)

Information not available

SECTION 8. Exposure controls/personal protection

8.1. Control parameters

**DIPOTASSIUM HEXAFLUORO-
TITANATE**

--	--	--	--	--	--	--	--	--



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 6/15

Predicted no-effect concentration - PNEC								
Normal value in fresh water	0,02	mg/l						
Normal value in marine water	0,002	mg/l						
Normal value for fresh water sediment	24,45	mg/kg/d						
Normal value for marine water sediment	4,89	mg/kg/d						
Normal value for marine water, intermittent release	0,11	mg/l						
Normal value of STP microorganisms	1,5	mg/l						
Normal value for the terrestrial compartment	19,1	mg/kg/d						
Health - Derived no-effect level - DNEL / DMEL								
	Effects on consumers				Effects on workers			
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic systemic	Acute local	Acute systemic	Chronic local	Chronic systemic
Inhalation						5,2 mg/m ³	5,2 mg/m ³	5,2 mg/m ³
Skin		37,5 mg/kg bw/d		37,5 mg/kg bw/d		75 mg/kg bw/d		75 mg/kg bw/d

DIATOMITE								
Predicted no-effect concentration - PNEC								
Normal value of STP microorganisms	100	mg/l						
Health - Derived no-effect level - DNEL / DMEL								
	Effects on consumers				Effects on workers			
Route of exposure	Acute local	Acute systemic	Chronic local	Chronic systemic	Acute local	Acute systemic	Chronic local	Chronic systemic



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 7/15

Oral				18,7 mg/ kg bw/d				
Inhalation				0,05 mg/ m3				0,05 mg/ m3

VND = hazard identified but no DNEL/PNEC available ; NEA = no exposure expected ; NPI = no hazard identified ; LOW = low hazard ; MED = medium hazard ; HIGH = high hazard.

During the risk assessment process, it is essential to take into consideration the ACGIH occupational exposure levels for inert particulate not otherwise classified (PNOC respirable fraction: 3 mg/m³; PNOC inhalable fraction: 10 mg/m³). For values above these limits, use a P type filter, whose class (1, 2 or 3) must be chosen according to the outcome of risk assessment.

8.2. Exposure controls

As the use of adequate technical equipment must always take priority over personal protective equipment, make sure that the workplace is well aired through effective local aspiration.

When choosing personal protective equipment, ask your chemical substance supplier for advice.

Personal protective equipment must be CE marked, showing that it complies with applicable standards.

Provide an emergency shower with face and eye wash station.

Exposure levels must be kept as low as possible to avoid significant build-up in the organism. Manage personal protective equipment so as to guarantee maximum protection (e.g. reduction in replacement times).

HAND PROTECTION

In the case of prolonged contact with the product, protect the hands with penetration-resistant work gloves (see standard EN 374).

Work glove material must be chosen according to the use process and the products that may form. Latex gloves may cause sensitivity reactions.

SKIN PROTECTION

Wear category II professional long-sleeved overalls and safety footwear (see Regulation 2016/425 and standard EN ISO 20344).

Wash body with soap and water after removing protective clothing.

EYE PROTECTION

Wear airtight protective goggles (see standard EN 166).

RESPIRATORY PROTECTION

Use a type P filtering facemask, whose class (1, 2 or 3) and effective need, must be defined according to the outcome of risk assessment (see standard EN 149).

ENVIRONMENTAL EXPOSURE CONTROLS

The emissions generated by manufacturing processes, including those generated by ventilation equipment, should be checked to ensure compliance with environmental standards.

SECTION 9. Physical and chemical properties

9.1. Information on basic physical and chemical properties

Properties	Value	Information
------------	-------	-------------



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 8/15

Appearance	powder	
Colour	blue	
Odour	perfumed	
Melting point / freezing point	not available	Reason for missing data: not relevant for the characterization of the product
Initial boiling point	not available	Reason for missing data: not relevant for the characterization of the product
Flammability	not flammable	
Lower explosive limit	not applicable	Reason for missing data: the substance/mixture is a solid
Upper explosive limit	not applicable	Reason for missing data: the substance/mixture is a solid
Flash point	not applicable	Reason for missing data: the substance/mixture is a solid
Auto-ignition temperature	not applicable	Reason for missing data: the substance/mixture is a solid
Decomposition temperature	not available	Reason for missing data: not relevant for the characterization of the product
pH	not applicable	Reason for missing data: substance/mixture is non-soluble (in water)
Kinematic viscosity	not applicable	Reason for missing data: the substance/mixture is a solid
Solubility	insoluble in water	
Partition coefficient: n-octanol/water	not applicable	Reason for missing data: substance/mixture is non-soluble (in water)
Vapour pressure	not available	Reason for missing data: not relevant for the characterization of the product
Density and/or relative density	2,3 g/cm ³	
Relative vapour density	not applicable	Reason for missing data: the substance/mixture is a solid
Particle characteristics	not available	Reason for missing data: not relevant for the characterization of the product

9.2. Other information

9.2.1. Information with regard to physical hazard classes

Information not available

9.2.2. Other safety characteristics

Information not available



Distributed by:
Omnident Dental- Handelsgesell-
schaft mbH
Gutenbergring 5, 63110 Rodgau



SECTION 10. Stability and reactivity

10.1. Reactivity

There are no particular risks of reaction with other substances in normal conditions of use.

10.2. Chemical stability

The product is stable in normal conditions of use and storage.

10.3. Possibility of hazardous reactions

The powders are potentially explosive when mixed with air.

10.4. Conditions to avoid

Avoid environmental dust build-up.

10.5. Incompatible materials

Information not available

10.6. Hazardous decomposition products

Information not available

SECTION 11. Toxicological information

In the absence of experimental data for the product itself, health hazards are evaluated according to the properties of the substances it contains, using the criteria specified in the applicable regulation for classification.

It is therefore necessary to take into account the concentration of the individual hazardous substances indicated in section 3, to evaluate the toxicological effects of exposure to the product.

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Metabolism, toxicokinetics, mechanism of action and other information

Information not available

Information on likely routes of exposure

Information not available

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Information not available

Interactive effects



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 10/15

Information not available

ACUTE TOXICITY

ATE (Inhalation) of the mixture:	Not classified (no significant component)
ATE (Oral) of the mixture:	>2000 mg/kg
ATE (Dermal) of the mixture:	Not classified (no significant component)

DIPOTASSIUM HEXAFLUOROTITANATE

LD50 (Oral):	324 mg/kg OECD 401
--------------	--------------------

DIATOMITE

LD50 (Oral):	> 2000 mg/kg OECD 401
--------------	-----------------------

SKIN CORROSION / IRRITATION

Does not meet the classification criteria for this hazard class

SERIOUS EYE DAMAGE / IRRITATION

Causes serious eye irritation

RESPIRATORY OR SKIN SENSITISATION

Does not meet the classification criteria for this hazard class

GERM CELL MUTAGENICITY

Does not meet the classification criteria for this hazard class

CARCINOGENICITY

Does not meet the classification criteria for this hazard class

REPRODUCTIVE TOXICITY

Does not meet the classification criteria for this hazard class

STOT - SINGLE EXPOSURE

Does not meet the classification criteria for this hazard class

STOT - REPEATED EXPOSURE

May cause damage to organs

ASPIRATION HAZARD

Does not meet the classification criteria for this hazard class

11.2. Information on other hazards



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 11/15

Based on the available data, the product does not contain substances listed in the main European lists of potential or suspected endocrine disruptors with human health effects under evaluation.

SECTION 12. Ecological information

Use this product according to good working practices. Avoid littering. Inform the competent authorities, should the product reach waterways or contaminate soil or vegetation.

12.1. Toxicity

DIPOTASSIUM HEXAFLUOROTITANATE		
LC50 - for Fish		172 mg/l/96h OECD 203
EC50 - for Crustacea		48,2 mg/l/48h OECD 203
EC50 - for Algae / Aquatic Plants		11 mg/l/72h OECD 211
EC10 for Algae / Aquatic Plants		1,4 mg/l/72h OECD 201
Chronic NOEC for Crustacea		10 mg/l OECD 211
DIATOMITE		
LC50 - for Fish		> 100 mg/l/96h OECD 203
EC50 - for Crustacea		> 100 mg/l/48h OECD 202
EC50 - for Algae / Aquatic Plants		> 100 mg/l/72h OECD 201

12.2. Persistence and degradability

DIPOTASSIUM HEXAFLUOROTITANATE		
Degradability: information not available		
DIATOMITE		
Degradability: information not available		

12.3. Bioaccumulative potential

Information not available



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 12/15

12.4. Mobility in soil

Information not available

12.5. Results of PBT and vPvB assessment

On the basis of available data, the product does not contain any PBT or vPvB in percentage $\geq$ than 0,1%.

12.6. Endocrine disrupting properties

Based on the available data, the product does not contain substances listed in the main European lists of potential or suspected endocrine disruptors with environmental effects under evaluation.

12.7. Other adverse effects

Information not available

SECTION 13. Disposal considerations

13.1. Waste treatment methods

Reuse, when possible. Product residues should be considered special hazardous waste. The hazard level of waste containing this product should be evaluated according to applicable regulations.

Disposal must be performed through an authorised waste management firm, in compliance with national and local regulations.

CONTAMINATED PACKAGING

Contaminated packaging must be recovered or disposed of in compliance with national waste management regulations.

SECTION 14. Transport information

The product is not dangerous under current provisions of the Code of International Carriage of Dangerous Goods by Road (ADR) and by Rail (RID), of the International Maritime Dangerous Goods Code (IMDG), and of the International Air Transport Association (IATA) regulations.

14.1. UN number or ID number

not applicable

14.2. UN proper shipping name

not applicable

14.3. Transport hazard class(es)

not applicable



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 13/15

14.4. Packing group

not applicable

14.5. Environmental hazards

not applicable

14.6. Special precautions for user

not applicable

14.7. Maritime transport in bulk according to IMO instruments

Information not relevant

SECTION 15. Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Seveso Category - Directive 2012/18/EU: None

Restrictions relating to the product or contained substances pursuant to Annex XVII to EC Regulation 1907/2006
None

Regulation (EU) 2019/1148 - on the marketing and use of explosives precursors
not applicable

Substances in Candidate List (Art. 59 REACH)
On the basis of available data, the product does not contain any SVHC in percentage $\geq$ than 0,1%.

Substances subject to authorisation (Annex XIV REACH)
None

Substances subject to exportation reporting pursuant to Regulation (EU) 649/2012:
None

Substances subject to the Rotterdam Convention:
None

Substances subject to the Stockholm Convention:
None

Healthcare controls



LASCOD S.p.A.
OMNI ALGINAT BLAU
PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnident Dental- Handelsgesellschaft mbH
Gutenbergring 5, 63110 Rodgau



Page 14/15

Workers exposed to this chemical agent must not undergo health checks, provided that available risk-assessment data prove that the risks related to the workers' health and safety are modest and that the 98/24/EC directive is respected.

15.2. Chemical safety assessment

A chemical safety assessment has not been performed for the preparation/for the substances indicated in section 3.

SECTION 16. Other information

Text of hazard (H) indications mentioned in section 2-3 of the sheet:

Acute Tox. 4	Acute toxicity, category 4	
STOT RE 2	Specific target organ toxicity - repeated exposure, category 2	
Eye Dam. 1	Serious eye damage, category 1	
Eye Irrit. 2	Eye irritation, category 2	
H302	Harmful if swallowed.	
H373	May cause damage to organs through prolonged or repeated exposure.	
H318	Causes serious eye damage.	
H319	Causes serious eye irritation.	

LEGEND:

- ADR: European Agreement concerning the carriage of Dangerous goods by Road
- ATE: Acute Toxicity Estimate
- CAS: Chemical Abstract Service Number
- CE50: Effective concentration (required to induce a 50% effect)
- CE: Identifier in ESIS (European archive of existing substances)
- CLP: Regulation (EC) 1272/2008
- DNEL: Derived No Effect Level
- EmS: Emergency Schedule
- GHS: Globally Harmonized System of classification and labeling of chemicals
- IATA DGR: International Air Transport Association Dangerous Goods Regulation
- IC50: Immobilization Concentration 50%
- IMDG: International Maritime Code for dangerous goods
- IMO: International Maritime Organization
- INDEX: Identifier in Annex VI of CLP
- LC50: Lethal Concentration 50%
- LD50: Lethal dose 50%
- OEL: Occupational Exposure Level
- PBT: Persistent bioaccumulative and toxic as REACH Regulation
- PEC: Predicted environmental Concentration
- PEL: Predicted exposure level
- PNEC: Predicted no effect concentration
- REACH: Regulation (EC) 1907/2006
- RID: Regulation concerning the international transport of dangerous goods by train
- TLV: Threshold Limit Value



LASCOD S.p.A.

OMNI ALGINAT BLAU PFEFFERMINZ

Revision nr. 1
Dated 19/09/2023
First compilation
Printed on 19/09/2023

Distributed by:
Omnicent Dental- Handelsgesell-
schaft mbH
Gutenbergring 5, 63110 Rodgau



Page 15/15

- TLV CEILING: Concentration that should not be exceeded during any time of occupational exposure.
- TWA: Time-weighted average exposure limit
- TWA STEL: Short-term exposure limit
- VOC: Volatile organic Compounds
- vPvB: Very Persistent and very Bioaccumulative as for REACH Regulation
- WGK: Water hazard classes (German).

GENERAL BIBLIOGRAPHY

1. Regulation (EC) 1907/2006 (REACH) of the European Parliament
 2. Regulation (EC) 1272/2008 (CLP) of the European Parliament
 3. Regulation (EU) 2020/878 (II Annex of REACH Regulation)
 4. Regulation (EC) 790/2009 (I Atp. CLP) of the European Parliament
 5. Regulation (EU) 286/2011 (II Atp. CLP) of the European Parliament
 6. Regulation (EU) 618/2012 (III Atp. CLP) of the European Parliament
 7. Regulation (EU) 487/2013 (IV Atp. CLP) of the European Parliament
 8. Regulation (EU) 944/2013 (V Atp. CLP) of the European Parliament
 9. Regulation (EU) 605/2014 (VI Atp. CLP) of the European Parliament
 10. Regulation (EU) 2015/1221 (VII Atp. CLP) of the European Parliament
 11. Regulation (EU) 2016/918 (VIII Atp. CLP) of the European Parliament
 12. Regulation (EU) 2016/1179 (IX Atp. CLP)
 13. Regulation (EU) 2017/776 (X Atp. CLP)
 14. Regulation (EU) 2018/669 (XI Atp. CLP)
 15. Regulation (EU) 2019/521 (XII Atp. CLP)
 16. Delegated Regulation (UE) 2018/1480 (XIII Atp. CLP)
 17. Regulation (EU) 2019/1148
 18. Delegated Regulation (UE) 2020/217 (XIV Atp. CLP)
 19. Delegated Regulation (UE) 2020/1182 (XV Atp. CLP)
 20. Delegated Regulation (UE) 2021/643 (XVI Atp. CLP)
 21. Delegated Regulation (UE) 2021/849 (XVII Atp. CLP)
 22. Delegated Regulation (UE) 2022/692 (XVIII Atp. CLP)
- The Merck Index. - 10th Edition
 - Handling Chemical Safety
 - INRS - Fiche Toxicologique (toxicological sheet)
 - Patty - Industrial Hygiene and Toxicology
 - N.I. Sax - Dangerous properties of Industrial Materials-7, 1989 Edition
 - IFA GESTIS website
 - ECHA website
 - Database of SDS models for chemicals - Ministry of Health and ISS (Istituto Superiore di Sanità) - Italy

Note for users:

The information contained in the present sheet are based on our own knowledge on the date of the last version. Users must verify the suitability and thoroughness of provided information according to each specific use of the product.

This document must not be regarded as a guarantee on any specific product property.

The use of this product is not subject to our direct control; therefore, users must, under their own responsibility, comply with the current health and safety laws and regulations. The producer is relieved from any liability arising from improper uses.

Provide appointed staff with adequate training on how to use chemical products.

CALCULATION METHODS FOR CLASSIFICATION

Chemical and physical hazards: Product classification derives from criteria established by the CLP Regulation, Annex I, Part 2. The data for evaluation of chemical-physical properties are reported in section 9.

Health hazards: Product classification is based on calculation methods as per Annex I of CLP, Part 3, unless determined otherwise in Section 11.

Environmental hazards: Product classification is based on calculation methods as per Annex I of CLP, Part 4, unless determined otherwise in Section 12.